

## 定性濾紙 Qualitative Filter Papers



· 定性分析，過濾結晶沈澱物

編號	No. 1	No. 2	No. 231	No. 232	No. 131	No. 235	No. 101	※ No. 4A
直徑尺寸(mm)	55,70,90,110,125,150,185,240,285,300,330,360,400,500,600,600x600							
包裝	100張/盒							

※No.4A為硬質濾紙

## 定量濾紙 Quantitative Filter Papers



· 重量分析，影響環境物質監測

編號	No. 3	No. 5A	No. 5B	No. 5C	No. 6	No.7
直徑尺寸(mm)	55,70,90,110,125,150,185,240,285,300,330,360,400,500,600					
包裝	100張/盒					

## 色層用分析濾紙 Chromatography Papers



· 色層分析，電泳及墨點分析

尺寸(mm)/編號	No. 50	No. 51A	No. 51B	No. 514A	No. 526	No. 590
20x400	100張/盒					
200x200	50張/盒					
600x600	50張/盒					

## 圓筒濾紙 Extraction Thimbles



· 高純度無縫筒狀濾紙，應用壓油、殺蟲劑、污染物分析、空氣粉塵採樣

編號	外直徑(mm) x 高(mm)	包裝
No. 84	20x90, 22x90, 25x90, 28x100, 30x100, 33x120, 35x120, 40x150, 45x150, 53x150, 60x200, 75x210	25支/盒
No. 86R		
No. 88R	20x90, 22x90, 28x100, 30x100, 33x120, 35x120, 40x150, 45x150, 53x150, 60x200, 75x210	10支/盒
No. 88RH	25x90	
No. 89		
No. 88R		
No. 89S		

## 玻璃纖維濾紙 Glass Fiber Filters



· 預過濾用，以延長濾膜使用時間 · 耐熱溫度：500℃  
· 水空氣污染分析 · 液體淨化

編號	DP-70	GA-55	GA-100	GA-200	GB-100R	GB-140	GC-50	GC-90	GD-120	GF-75	GS-25
包裝	50張/盒	100張/盒	100張/盒	50張/盒	100張/盒	100張/盒	100張/盒	100張/盒	50張/盒	100張/盒	100張/盒
直徑尺寸(mm)	21, 24, 25, 26, 37, 45, 47, 55, 70, 90, 110, 125, 150										

※另有方型濾紙300x300mm，包裝：10張/盒

## 矽纖維濾紙 Quartz Fiber Filters



· 高溫(>1000℃) 酸氣取樣、空氣污染分析

編號	包裝	直徑尺寸(mm)
QR-100	100張/盒	21, 24, 26, 37, 45, 47, 55, 70, 90, 110, 125, 150
QR-200	50張/盒	

● 另有 No.545 新生兒採血濾紙 PKU Filter Papers, 350mmx540mm, 包裝：500張/盒