

# Laboratory Shaker (Horizontal Shaking)

MK201D

## Rotary and reciprocate shaking motion



- Compact, space-saving design
- Changeable rotary and reciprocate motion for mixing, extracting and stirring of samples
- Digital display of shaking frequency
- Equipped with shaking timer
- Easy to assemble and remove accessories
- Various shaking modes when used with different shaking stages and racks
- Can be placed inside IN-602CSW incubator for shaking incubation

### Specifications

Model	MK201D
Shaking mode	Rotary and Reciprocate (Manual operation)
Shaking range	Rotary: 30mm Reciprocate: 30mm
Shaking frequency	20~200rpm, Random adjustment
Display type	Digital
Timer	0.5sec.~100hrs.
External dimensions	W442×D415×H130mm
Weight	19kg
Power source (50/60Hz)	AC115V Approx. 2A / AC220V Approx. 1.3A

### Operational accessories

Shaking rack for erlenmeyer flask clamp, Fixed type



Product code	Capacity	No. of clamps
232170	100ml	20pcs.
232171	200ml	10pcs.
232172	300ml	6pcs.
232173	500ml	4pcs.
232174	1,000ml	2pcs.

Shaking rack for erlenmeyer flask clamp, Inclined type



2 sets example (Inclination stage sold separately)

For 1 set		
Product code	Capacity	No. of clamps
232175	100ml	8pcs.
232176	200ml	4pcs.
232177	500ml	2pcs.

Shaking rack for spitz tube



2 sets example (Inclination stage sold separately)

For 1 set	
Size 15ml	12pcs.
E.D. ø16.5 x L110mm or less	
Product code: 232178	

Shaking rack for centrifugal tube



2 sets example (Inclination stage sold separately)

For 1 set	
Size 50mL	8pcs.
E.D. ø30×L110mm or less	
Product code: 232179	

Inclination setting stage



Necessary in tilt cases, up to 2 sets  
Product code: 232078

Shaking rack with sticky non-skid sheet



W400×D330mm Material: silicon rubber  
Product code: 232182

(For flat bottom flask use only)

Shaking rack for SUS spring almighty type A




A type: for polyethylene bottle
Size 250ml up to 9pcs.
Product code: 232180

Shaking rack for SUS spring almighty type B



B type: for erlenmeyer flask and test tube
17 x 15 = 255 grids
Product code: 232181

Two layer shaking rack for SUS spring almighty type C



C type: for erlenmeyer flask and test tube
For erlenmeyer flask
100ml×16pcs, 200ml×10pcs,
500ml×4pcs, 1,000ml×2pcs.
For test tube ø16mm×105pcs. when 45°inclination
Product code: 231398

\*Glassware and polyethylene bottles not included.